

# BILGATH KIT

PERCUTANEOUS DRAINAGE CATHETER WITH PIGTAIL POINT WITH OR WITHOUT LOCKING SYSTEM

## LINEA BILGATH

Complete kit formed by:

- Polyethylene Bilcath Pigtail catheter
- n° 1 Ecojekt BS needle
- n° 1 Ecojekt needle
- Straight nitinol guidewire .038
- Dilators
- Fixation disc
- Bag connection tube
- Syringe connection tube
- Mini scalpel
- 3way stopcock
- Drape
- Sterile towel
- Washing valve
- Connection tube

without locking system

CODICE	Diametro	Lunghezza	Box
BILK0850-PN	Fr.8	500mm	1
BILK1050-PN	Fr.10	500mm	1
BILK1250-PN	Fr.12	500mm	1

with locking system


CODICE	Diametro	Lunghezza	Box
BILK0850-PL	Fr.8	500mm	1
BILK1050-PL	Fr.10	500mm	1
BILK1250-PL	Fr.12	500mm	1

Legenda:

P = Pigtail

N = without locking system

L = with locking system



PERCUTANEOUS DRAINAGES