

INTRODUCER NEEDLE FOR COAXIAL BIOPSY "N" TYPE

INTRODUCER NEEDLE FOR PERCUTANEOUS PROCEDURES

H.S. Hospital Service S.p.A. è una società certificata ISO 9001 e EN 46001

Introducer needle for percutaneous procedures

Features:

- *Blunt and tapered pointed needle*
- *Echogenic marker*
- *Depth marks*
- *Clear luer hub*
- *Stylet with colour code and pyramidal point*

CODE	Diameter	Box
IN1005N	10G x 50mm	10
IN10075N	10G x 75mm	10
IN1305N	13G x 50mm	10
IN13075N	13G x 75mm	10
IN13095N	13G x 95mm	10
IN13110N	13G x 110mm	10
IN13140N	13G x 140mm	10
IN1505N	15G x 50mm	10
IN15075N	15G x 75mm	10
IN15095N	15G x 95mm	10
IN15110N	15G x 110mm	10
IN15140N	15G x 140mm	10
IN1705N	17G x 50mm	10
IN17075N	17G x 75mm	10
IN17095N	17G x 95mm	10
IN17140N	17G x 140mm	10



CE 0373

HISTOLOGICAL BIOPSY



H.S. HOSPITAL SERVICE S.p.A.

Via delle Valli snc - 04011 Aprilia (LT)
Tel. 069201961 - Fax 0692727871 - Fax Uff. ordini 069275519
www.hshospitalservice.com e-mail: hs@hshospitalservice.com