

PRECISA VACUUM SET

PRECISA VACUUM SYSTEM

- Device for vacuum creation - Trolley included
- 4 Shelves
- Single pedal

CODE
MVS-S

PRECISA VACUUM SET

- Tru-cut point needle
- Protective cap
- ABS Syringe
- Connection Tube
- Removable tru-cut stylet with whistle point and 20mm notch
- Luer Lock Hub
- A two rings handle which allows to perform the biopsy using only one hand



CODE SIZES

MVSI110	11G x 100mm with introducer needle
MVSI115	11G x 150mm with introducer needle
MVSI410	14G x 100mm
MVSI415	14G x 150mm
MVSI420	14G x 200mm
MVSI610	16G x 100mm
MVSI615	16G x 150mm

Work in progress



HS HOSPITAL SERVICE S.p.A.
Via delle Valli snc - 04011 Aprilia (LT)
Tel. +39 06 9201961-Fax +39 0692727871
www.hshospitalservice.com